


Quick Start Guide

Prolaio Bracelet



This guide provides step-by-step instructions
to use the Prolaio system.

prolaio


Your kit contains the following supplies:

- A. Corsano bracelet
- B. Bracelet charger and cable



Connect your bracelet


- Tap the screen. *The Home Page appears with Action Cards to guide you with next steps.*
- Tap the **BRACELET** tab next to **PATCH** tab on the Status card.
- Tap **Open Connection Guide** and follow the directions in the app.



STEP 1

Power-On the Bracelet

- Please connect the charger provided in your monitoring kit to the bracelet. **It's important to ensure that it is plugged into an outlet.**
- The LEDs on the side of the device will start blinking, indicating that it is ready to be connected.
- Tap **Next** on the Prolaio phone to proceed to the next step.




STEP 2

Pairing the Bracelet

- Wait a few seconds while the app searches for your bracelet.
- Once found, a dialog will prompt you to remove the bracelet from the charger and wait for the app to establish the connection.
- Tap **Pair** on the Bluetooth pairing request dialogue.
- The app automatically updates the bracelet's software whenever necessary. A percentage of update progress is presented.

Please remove the bracelet from the charger



58
sec

Please remove the charger and wait for the connection to be established

STEP 3

Wrist selection

- On the app, choose the wrist where you will wear the bracelet.
- This selection will help you position the bracelet LED orientation.
- Tap **Next** on the Prolaio phone to proceed to the next step.


Step 3 of 5

Select wrist for your bracelet.

LEFT HAND

RIGHT HAND


The LED orientation on the bracelet should face away from you when worn on the **left** hand.



STEP 4

Bracelet position

- The bracelet should be placed about 1 inch above the bone of your wrist, as illustrated below.
- The LED orientation on the bracelet should face towards you when worn on the right hand and away from you when worn on the left hand.
- Please adjust the strap of your bracelet so that it fits snugly but not too tightly. This will help you keep the bracelet in place when you use your hands without compromising comfort. Please ensure that the light emitted from the sensor is not visible between your skin and the bracelet.
- Tap **Next** on the Prolaio phone to proceed to the next step.



STEP 5

Final confirmation

- Finally, the app will confirm that the bracelet is worn correctly and the app is receiving its data.

Step 5 of 5




Confirming that bracelet is being worn.

CAUTION

Do not use the Bracelet device during an MRI scan or in a location where it will be exposed to strong electromagnetic forces.

If discomfort or irritation occurs, the Bracelet should be removed. If mild soreness or redness is experienced after removing the Bracelet, Do not wear the bracelet on the same wrist.

Icon and Description

	Manufacturer		Humidity Limitation
	Use-by date		Caution
	Lot Number		Prescription Use Only
	Temperature Limit		Consult instructions for use Operator's manual
	MR Unsafe		

prolaio

Contact us

(855) 869-9054

Monday-Friday, 9am - 5pm CST
(excluding major holidays)

support@prolaio.com

230 W. Monroe, Suite 2560
Chicago, IL 60606 USA



Guidelines for best use

- You can wear the bracelet when showering, but it is not advisable due to possible skin irritation.
- Do not cover, submerge it in water, or use it in a sauna.
- The bracelet's battery can last for up to 3½ days. After three days, the phone will alert you when the battery is low, letting you know that you will need to charge the bracelet soon.
- It can take up to 2 hours to fully charge the bracelet. Therefore, recharging the bracelet two hours before bedtime on the third day of wearing it is recommended.
- The Prolaio phone should be kept close to and in the same room as you as much as possible. 20 feet (6 meters) or less is ideal.
- The Prolaio phone and bracelet can be out of range of one another for up to 15 hours (MAX) per day without losing data.
- The Prolaio phone should always be kept powered on and charged.
- You can travel with the bracelet. When passing through airport security, remember to take it off your wrist and put it back on after completing the security check. Additionally, remember to turn off the phone when flying.
- If you plan to leave the country, inform the professional who issued the Prolaio phone to you.
- All components and accessories are magnetic resonance (MR) unsafe and can pose a projectile hazard in the MR environment, and therefore, must be kept out of the magnetic resonance imaging (MRI) scanner room.
- **DO NOT APPLY THE DEVICE TO ANY REDDENED, IRRITATED, OR OPEN LESION ON YOUR SKIN. NOTIFY YOUR HEALTHCARE PROFESSIONAL IF ANY SKIN IRRITATION OCCURS.**